Gebrauchsinformation

DE

Pulp Protection

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Technical Data ProHALO

Complies with ISO 4049

Range of dimensions of organic filler particles:

0.04 - 2.9 µm

Inorganic filler content by volume:

56 %

Inorganic filler content by weight:

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• Lining and repair of ceramic and composite restorations

• More translucent material in the layer technique

Contra-indications

If the site cannot be isolated after enamel etching and during application and curing of ProHALO. If allergic to components of ProHALO. If oral hygiene is poor.

Side effects

No systemic side effects are known. Contact allergies with products of similar composition have been reported in isolated cases.

Interactions with other agents

Agents which contain eugenol and/or clove oil may affect the polymerization of ProHALO. The use of zinc-oxide-eugenol cements in combination with ProHALO should therefore be avoided. Discolorations may occur when using cationic mouth rinses, plaque indicators or chlorhexidine.

Application

Cleansing and isolation

Cleanse the tooth to treat and the adjacent teeth.

Select color and create a dry working area.

Preparing the cavity

When preparing the cavity, every effort should be made to use a technique that is gentle to the tooth substance (principle of adhesive Restoration Technique). Enamel edges in the anterior region should be beveled. In the posterior region, only sharp enamel edges should be slightly beveled or rounded.

Pulp Protection

A suitable pulp-protective liner should be applied in close proximity to the pulp.

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IT

Pulp Protection

A suitable pulp-protective liner should be applied in close proximity to the pulp.

Procedure according to manufacturer's instructions for use.

Application of bonding agent

Apply the bonding agent. Procedure according to manufacturer's instructions for use.

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Notas

For the treatment of caries, dentine, dental cement and composite restorations.

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