

## **Data Sheet ROEKO Cotton rolls**

General deviceThe dental rolls are cylindrically shaped cotton wool products which are available in<br/>various sizes and are suited for use with adults and children.

The dental rolls are supplied unsterile. Due to the production process they are low of germs and they are intended for single use.

The dental rolls are brought inside the mouth for a brief time and come into contact with lips, oral mucosa and tooth surfaces.

Because the dental rolls are uniformly manufactured to length, they can be hygienically stored in compatible dispenser systems.

In these roll dispensers, the dental rolls are removed on a first in, first out principle.

Intended use: The dental rolls are intended to be used during dental treatments that require moisture control/tooth isolation. They absorb moisture and retract cheek, lips and tongue to keep operative areas accessible and free of excess moisture.

The dental rolls are available in the following sizes:

Art. No.	Size	Ø/mm	Length/mm
60024049	1	8,0 ±0,3	37±1,0
60024050	2	10,0 ±0,3	37±1,0
60024051	3	12,0 ±0,3	37±1,0

Manufacturer: Dent4you GmbH

Storage and handling:intended for single use only,<br/>medical device,<br/>keep dry,<br/>do not use after the expiration date has passed

**Composition:** cotton fibers, cellulose, cellulose wadding 2-ply oxygen bleached, adhesive

Shelf life: 5 years