

## Data Sheet Guttapercha Points

**General device description:** Guttapercha Points are used in combination with sealers to obturate the root canal after root canal treatment.

The Guttapercha Points are available in different ISO sizes, in different tapers and with a top color coding, full color coding or in pink. Also conventional sizes are available to offer the dentist more options for various condensation techniques.

**Intended use:** Guttapercha Points are used as a filling material (permanent use) of the cleaned tooth root canals to fill the root canal of a prepared tooth, in strict accordance with internationally recognized standards. In addition, it is only intended for use by qualified medical personnel trained in its use.

The Guttapercha Points are available in the following sizes:

Art. No	Description	Size
60024220	Gutta percha Prodonto Points	15
60024221	Gutta percha Prodonto Points	20
60024224	Gutta percha Prodonto Points	25
60024226	Gutta percha Prodonto Points	30
60024227	Gutta percha Prodonto Points	35
60024228	Gutta percha Prodonto Points	40
60024229	Gutta percha Prodonto Points	45
60024230	Gutta percha Prodonto Points	50
60024231	Gutta percha Prodonto Points	55
60024232	Gutta percha Prodonto Points	60
60024234	Gutta percha Prodonto Points assorted	15-40
60024235	Gutta percha Prodonto Points assorted	45-80

**Manufacturer:** Dent4you GmbH

**ISO standard:** according to ISO 6887

---

**Storage and handling:** medical device,  
intended for single use only,  
keep dry,  
do not use after the expiration date has passed,  
protect from sunlight,  
storage below 25°C and 77°F

**Composition:** Guttapercha mixture, colouring agents

**Shelf life:** 5 years