

Data Sheet Paper Points

General device description: Paper Points are made from highly absorbent paper intended to clean the root canal of a tooth by insertion and removal supplied in sterilized condition.

Intended use: Paper Points are intended to dry the root canal after preparation and irrigation.

The Paper Points are available in the following sizes:

Art. No	Description	Size
60024236	Paper Points Prodonto	15
60024237	Paper Points Prodonto	20
60024238	Paper Points Prodonto	25
60024239	Paper Points Prodonto	30
60024240	Paper Points Prodonto	35
60024241	Paper Points Prodonto	40
60024242	Paper Points Prodonto	45
60024243	Paper Points Prodonto	50
60024244	Paper Points Prodonto	55
60024245	Paper Points Prodonto	60
60024246	Paper Points Prodonto	70
60024247	Paper Points Prodonto	80
60024249	Paper Points Prodonto assorted	15-40
60024250	Paper Points Prodonto assorted	45-80

Manufacturer: Dent4you GmbH

Storage and handling: medical device,
intended for single use only,
keep dry,
do not use after the expiration date has passed,
do not use if packaging is damaged

Composition: Paper, Water, Color

Shelf life: 5 years